

ONE SHADING MEETS BULK FILLING



NEW



Charisma Bulk Flow ONE.

The first flowable one-shade bulk-filling composite for everyday cases.

Always grab the right shade: With the incredible shade-matching properties of Charisma Bulk Flow ONE the restoration becomes invisible by blending into the surrounding dentition.

Efficient handling: Cavities are filled quickly and easily with increments of up to 4mm without the need for an additional capping layer.

Long-lasting restorations: Low shrinkage stress, high flexural strength and wear resistance ensure durable restorations that withstand the mastication load in the posterior region.



Try it out for yourself: kulzer.com/try-charisma-bulk-flow-one



KULZER
MITSUI CHEMICALS GROUP

ONE LAYER FOR ALL SHADES



Charisma Bulk Flow ONE simplifies and streamlines the restorative process.

Charisma Bulk Flow ONE combines the self-levelling effect of a flowable material with a high curing depth and good stability after application. This allows to place increments of up to 4 mm inside a cavity.

With Charisma Bulk Flow ONE dentists can reduce treatment time without the need for shade selection or the placement of a universal composite as an additional capping layer.

With its range of indications from bulk filling through cavity lining to fillings in primary teeth, the efficient restoration procedure is further enhanced by the equally reliable and versatile bonding solution iBOND Universal.

VERSATILITY IN APPLICATIONS



EFFICIENT HANDLING



Excavated and prepared cavity under rubberdam-isolation.



Clinical follow-up after 4 weeks.



X-ray control: Very good radiopacity of Charisma Bulk Flow ONE

Art. code	Package	Content
Charisma Bulk Flow ONE – Refills		
66094725	CHARISMA BULK FLOW ONE SYR 1x2G	1 x 2g syringe 5 x cannulas 1 x IFU
Charisma Bulk Flow ONE – Value Packs		
66095336	Charisma BULK FLOW ONE VALUE-KIT SYRINGE	10 x 2g syringes 10 x 5 cannulas 1 x pictorial card 1 x IFU



Discover our new ONE: kulzer.com/try-charisma-bulk-flow-one